




AC- 673170	<b>CUSTOMER DRAWING</b>	Page 1 of 1
Revision 02	PS Microbatch IMP@CT™ Plate, 96 Well	 <b>greiner bio-one</b>
	Item-No. 673 170	



NOT TO SCALE



Cat.-No.	Description	A	B	C	D	D'	E	F	G	H	I
673 170	IMP@CT™ Plate	127.76 ±0.1	85.48 ±0.1	14.4 ±0.1	Ø 3.3	Ø 1.34	2.5	1.8	14.38	11.24	9.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date October 02, 2009	Date October 02, 2009	Date October 05, 2009	
Date November 7, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	