

EASYseal™, BREATHseal™, POWERseal™, VIEWseal™, SILVERseal™

Greiner Bio-One offers a full set of plate sealers for a wide range of routine and specialized applications, e.g. high-throughput screening, immunology, microbiology, molecular biology, and PCR.

EASYseal™ is recommended for all laboratory purposes, e.g. ELISA applications to protect samples from evaporation or contamination.

BREATHseal™ is a unique gas permeable sealer with a wide distribution of pores for excellent gas exchange. Available sterile and non sterile, this gas permeable sealer is suitable for a variety of applications including the culture of prokaryotic and eukaryotic cells.

POWERseal™ is coated with a pressure-sensitive silicon adhesive that will only adhere where pressure is applied. The textured film can be easily pierced with metal probes or pin tools.

VIEWseal™ is a special highly transparent sealer for precise optical measurements. The pressure-sensitive silicon coating will only adhere where pressure is applied.

SILVERseal™ is an extremely robust, heat resistant and pierceable aluminum film, coated with a pressure-sensitive acrylate adhesive.



Features

- Non-pyrogenic
- DNase- and RNase-free
- Free of DNA-polymerase-inhibitors
- Pierceable films
- Heat proofed up to +110°C
- Cold proofed down to -80°C

Ordering Information

Cat.-No.	Description	Qty/case
676 001	EASYseal™, adhesive sealer, standard, transparent	100
676 050	BREATHseal™, gas permeable, pierceable	50
676 051	BREATHseal™, gas permeable, pierceable, sterile	50
676 080	POWERseal™, advanced adhesive sealer, textured, pierceable	100
676 070	VIEWseal™, advanced adhesive sealer, transparent	100
676 090	SILVERseal™, adhesive sealer, aluminum, pierceable	100

Revision: August, 2003 – 073 065

www.gbo.com/bioscience

Germany (Main office)
Greiner Bio-One GmbH
Phone: (+49) 7022 948-0
E-Mail: info@de.gbo.com

Belgium
Greiner Bio-One N. V.
Phone: (+32) 2-4 61 09 10
E-Mail: info@be.gbo.com

France
Greiner Bio-One SAS
Phone: (+33) 169 -86 25 50
E-Mail: infos@fr.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Phone: (+81) 3-35 05-88 75
E-Mail: info@jp.gbo.com

Austria
Greiner Bio-One GmbH
Phone: (+43) 7583 6791-0
E-Mail: office@at.gbo.com

Netherlands
Greiner Bio-One B. V.
Phone: (+31) 172-42 09 00
E-Mail: info@nl.gbo.com

UK
Greiner Bio-One Ltd.
Phone: (+44) 1453-82 52 55
E-Mail: info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Phone: (+1) 704-261-78 00
E-Mail: info@us.gbo.com



greiner bio-one