

Cryo Accessories

A large choice of accessories for manipulation, storage and transportation of frozen samples



Puck for SPINE format sample holders

These pucks can be loaded with 10 frozen samples mounted on SPINE format holders.

Puck Supports, for transfer of SPINE sample holders

A precious tool for insertion and retrieval of the sample holders to/from the puck.



CX 100 Dry Shipper and its shipping case

The recommended dry shipper for carrying SPINE format pucks and other format.

Puck holder for storage Dewar

This pucks holder can accommodate five SPINE format pucks. Each puck can be placed and removed independently, in a very convenient way.

(For more details, see Bodin et al, J. Sync. Radiat. 2012 Mar;19(Pt2):288-9)



And many more...