



Simplify Storage & Shipping

Cryo-EM Pucks

Cryogenic Grid Box Storage and Shipping



- 12 wells per puck for Cryo Grid Boxes
- Unique identifier for each puck
- Indexed wells for sample tracking
- Tweezer slots for easy and secure box handling
- Accommodates pin type lid style boxes
- Shelved storage cane holds eight pucks
- Shelved shipping cane holds seven pucks
- Free custom puck engraving

➤ Easily Stored & Shipped



Simplify how you store and transport your cryogenically cooled Cryo-EM specimen grids. Each puck holds up to 24 round Cryo Grid Boxes (96 grids). An 8 puck capacity storage cane is compatible with standard refrigerators, such as the HC34. The compatible shelved shipping cane holds 7 pucks and can be used with a CX100 dryshipper.

mitegen.com/cem/empucks